



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG 15-EN MELATRONIC Autoclave (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1500 W,
- Has 4 working programs:
 - Universal program - 5.5 min sterilization time at 134 °C,
 - Quick program - sterilization time of 3.5 min at 134 °C,
 - Protective program - sterilization time 20.5 min at 121 °C,
 - Prion program - sterilization time 20.5 min at 134 °C,
- Capacity of sterilization chamber: 6.7 liters,
- Manual water infusion,
- Pressure release button,
- Menu in German,
- Included:
 - Power cord,
 - Tray stand,
 - Tray holder,
 - 2 pcs. - 11.5 x 34 cm trays,
- Dimensions: 44 x 50 x 33 cm,
- Weight: 19 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

