



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC LIFEPAK 20 Defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: through 1 built-in battery,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- Selectable language menu: PL, EN, FR, DE, ES, SV, IT, DA, PT, ZH, CS, RU, JA, HU, KO, NL, NO, FI,
- Lifepak 20 is a biphasic portable defibrillator, providing a complete response system for acute cardiac conditions,
- Technical specifications:
 - AED mode,
 - Synchronization,
 - ECG parameters monitoring,
 - Easy-to-read 6.2" display,
 - Built-in printer (paper width: 50 mm),
- Includes:
 - User manual in English (PDF),
 - Power cord,
 - ECG 10-lead cable,
 - AED cable,
- Dimensions: 26 x 24 x 21 cm,
- Weight: 6.5 kg,
- Cat D135,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

