



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



ATMOS C161 Medical suction unit (Ex-Display)

- New (ex-display),
- Real photos of the product,
- Made in Germany, (2020),
- Power supply: 230 V,
- Emergency power supply: through a new built-in 8.4 V / 6 Ah battery,
- Frequency: 50/60 Hz,
- Power: 45 VA,
- SSAK ATMOS C161 functional and comfortable medical suction machine for home and hospital use,
- The C 161 model of the suction unit is powered by the mains and the built-in battery, making it ideal for stationary care and the field. Its design and small size allow it to be used comfortably in small rooms or stored on shelves with small space,
- Technical specifications:
 - Max flow rate (L/min): 16/2,
 - Vacuum force adjustment: Infinitely adjustable to -80 kPa,
 - Serres tank with a capacity of 1 l.,
- Features:
 - Easy operation,
 - Maximum safety with triple protection against flooding of the device,
 - Low operating costs,
 - High efficiency,
 - Quiet operation,
 - Moisture monitoring system,
 - Control panel cover,
 - Built-in handle for easy transportation,
- Included:
 - User manual in English (PDF),
 - Power supply,
 - SERRES reservoir,
 - Suction-controlled suction tube with a length of 180 cm,
- Overall dimensions: 25.5 x 18 x 25 cm,
- Weight: 3.5 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

