



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VISION MEDICAL SMARTGRAFT Hair transplantation device (Ex-Display)

- New (ex-display),
- Realistic photos of the product,
- Made in the USA, (2019),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 770 VA,
- SmartGraft is a closed system for hair transplantation procedures. It allows hair collection, storage and hydration before reimplantation. The procedure is performed without a scalpel or the use of sutures,
- A highly efficient pump and a wide set of treatment heads allow hair extraction without damaging the hair. Hair survival rate reaches 98%. SmartGraft is equipped with a system for counting the extracted hairs,
- SmartGraft is one of the few hair thickening methods that provides such a high level of safety. During the procedure, the doctor has at his disposal 2 types of treatment needles of different diameters. Each is equipped with a special separator that does not allow the needle to enter the skin deeper than 1 millimeter. SmartGraft does not require the use of a scalpel, so no longitudinal scars form on the skin. It is a minimally invasive procedure performed in the office, rather than a complicated procedure,
- Indications:
 - Androgenetic alopecia in men and women,
 - Post-traumatic alopecia,
 - Baldness after burns,
 - Baldness after irradiation,
 - Baldness caused by inflammatory skin conditions,
 - Post-surgical scars,
 - Complementary to reconstructive scalp treatments,
 - Hair loss in the area of eyebrows, eyelashes, mustache and beard,
 - Reconstruction of chest hair,
- Includes:
 - SmartGraft Cart (cart),
 - Head,
 - Cable for the head,
 - PTFE filter - 2 pcs,
 - Footswitch - 2 pcs.
 - Needles - 5 pcs,
 - Tank set - 5 pcs,
 - Sprayer - 1 pcs,
 - Metal container - 3 pcs,
 - Set of screws for mounting the cart ,
 - Head holder (black),
 - Accessory holder (on shoulder),
 - Set of screws for accessory holder,
 - Head tip,
 - Power cord,
- Weight: 65.6 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

