



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ARJO ENCORE SARA PLUS Patient lift (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Sweden,
- Power supply: 230 V,
- Emergency power supply: through 4 new built-in 6 V / 5 Ah batteries,
- Frequency: 50/60 Hz,
- Sara Plus offers advanced options to support patient mobility more than any other instrument of its kind,
- Features:
 - Designed for operation by a single caregiver, Sara Plus can be a cost-effective solution for lifting, stepping and walking training,
 - PRO-ACTIVE - Active knee pad provides support and stability while allowing natural forward knee movement and ankle flexion during standing and sitting activities,
 - Detachable footrest - Detachable footrest allows stepping and walking practice with the patient's feet directly on the floor,
- Technical specifications:
 - Electric lift adjustment in the range: 90 - 150 cm,
 - Electric leg spacing adjustment: 64 - 110 cm,
- Includes:
 - User manual in English (PDF),
 - Power cord,
 - Power supply,
- Overall dimensions: 70 x 120 x 95 cm,
- Maximum load capacity: 190 KG,
- Weight: 79 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

