



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **EPENDORF CENTRIFUGE 5702 Laboratory centrifuge (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 200 W,
- The Centrifuge 5702 is a low-speed laboratory centrifuge specifically designed for use in cell culture and clinical research laboratories working with small to medium sample volumes,
- Designed for standard applications,
- Built-in backlit LCD display,
- Control panel keypad in English,
- Product features:
  - Compact design so the unit fits on any laboratory table,
  - Very quiet operation,
  - Low unit height for easy insertion and removal of samples,
  - Optimized for cell separation by gradient centrifugation,
  - "At set RPM" function starts the timer when the selected RPM is reached; to ensure repeatable centrifugation,
  - Push-button lock prevents unintentional parameter adjustment,
  - Stainless steel rotor chamber is easy to clean,
  - Electronic unbalance recognition for maximum safety,
  - Automatic lid lock,
  - Adjustable time: 0.5 - 99 minutes and continuous centrifugation,
- Technical specifications:
  - Max RCF: 3 000 × g,
  - Max speed: 4400 rpm,
  - Max capacity: 4 x 100 ml,
  - Noise level:
- Includes:
  - Rotor A-4-38,
  - 4 pcs. - 100 ml round basket,
  - 4 pcs. - adapter (for 2 15 ml conical tubes - 17.2 mm x 121 mm),
  - Instruction manual in Polish (PDF),
  - Power cord,
- Overall dimensions: 32 x 40 x 24/53 (with lid open) cm,
- Weight: 20 KG,
- Cat 05,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

