



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**DRAGER FABIUS TIRO Anesthesia machine (Respirator / Ventilator) (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 3.2 Ah batteries,
- Frequency: 50/60 Hz,
- Software: 3.37b,
- The Fabius Tiro is an inhalation anesthesia machine used in operating rooms, preoperative and recovery rooms,
- It has one port for vaporizer with DRAGER DW-2000 mounting system, no vaporizer included - can be retrofitted at additional cost,
- Powered by gases:
  - O2 connector,
  - N2O,
  - Air,
- Respirator driven by an electric motor (electronically controlled),
- Ventilator measurement functions:
  - P<sub>MAX</sub> (V) - 15-70 mbar,
  - V<sub>T</sub> - 20-1400 ml,
  - f - 4- 60 bpm,
  - Ti/Te - 4:1 to 1:4,
  - Tip/Ti - 0% - 50%,
  - PEEP - 0-20 mbar,
  - P<sub>insp</sub> (P) - 5-65 cmH2O,
  - InspFlow - 10-75 L/min,
  - PPS (SIMV/PS) - pressure support , from 3 to 20 cmH2O,
- Built-in LED monitor with O2, AIR, N2O gas measurement function, and control panel and ventilator function display,
- Easy to use and easy to read display,
- Infinitely adjustable gas flow rate - O2, AIR, N2O, from 0.0 to 12L/min, (digital displays + pressure gauges),
- Compact heated patient breathing system, located on a mobile arm, equipped with a dome and oxygen sensor,
- Built-in alarms, and low and high pressure valves,
- Control panel in Polish language,
- Selectable menus in the following languages: PL, BG, RO, HU, HR, CS, SK, RU, EN, FR, DE,
- Machine equipped with:
  - DRAGER INFINITY DELTA XL Cardiac monitor / patient monitor:
    - Reconditioned (used),
    - Technical condition: very good,
    - Visual condition: very good,
    - Real photos of the product,
    - Made in the USA,
    - Power supply: 230 V,
    - Emergency power supply via 1 new built-in 12 V / 2.3 Ah battery,
    - Frequency: 50/60 Hz,
    - Software: DELTA VF8.4-E,
    - Power: 70 W,
    - DRAGER INFINITY DELTA XL designed for multi-parameter patient monitoring,
    - Custom display options allow the operator to view and analyze data in graphical and numerical formats, compile real-time measurements and trends, and freely configure each screen element from curves to parameter names,
  - Equipped with:
    - 12.2" LCD color liquid crystal display
    - Navigation dial for easy information entry and screen navigation,
    - Device menu in Polish language,
    - Foiled control panel in English,
    - Convenient handle with and rugged housing for easy portability,
    - Pick and Go technologies,
    - Built-in tools such as hemodynamic calculators, drug calculators, screen cursors for dynamic curve analysis and more,
    - Audible and visual alarms,
  - Selectable menus in PL, ENG, CS, HUN,
  - DRAGER VAMOS PLUS gas module:
    - Reconditioned (used),
    - Technical condition: very good,
    - Visual condition: very good,
    - Real photos of the product,
    - Made in Germany,
    - Software: VAMOS 3.23,
    - DRAGER VAMOS PLUS gas module - used to reliably provide information on the concentration of carbon dioxide, nitrous oxide and inhaled anesthetics, in the patient's respiratory gas mixture,
    - Carbon dioxide concentration is shown on the screen in the form of a curve updated in real time,
    - Intuitive operation and easy-to-read screen,
    - Keypad on the control panel in English,
    - Comes with WaterLock Drager Wal water trap, providing greater protection against infection,
  - Selectable menu in GER, ENG, FR, IT, ESP, NL, SK,
  - DRAGER connection docking station:
    - Provides monitor mounting on stand, power supply, connection of additional screen and recorder, and data transfer,
- Included:
  - FABIUS TIRO user manual in Polish (PDF),
  - User manual for the DRAGER INFINITY DELTA XL patient monitor in Polish (PDF),
  - Brand new gas supply hoses with customer's choice of ends (AGA or DIN),
  - Brand new patient breathing system,
  - Brand new gas sampling line,
  - Power cord,
  - Wired modules:
    - MultiMed cable,
    - ECG 5-lead cable,
    - SPO2 sensor,
    - NiBP cuff (not included - can be retrofitted at additional cost),
    - Air tubing,
- Dimensions:
  - Overall: 90 x 86 x 150 cm,
  - 3 pull-out main drawers with dimensions:
    - 36.5 x 32.5 x 12.5 cm,
  - Pull-out shelf: 39 x 31.5 cm,
- Weight: 135 kg,
- Cat 09,
- In the domestic market (Poland), in the case of installation of the device within 30 days from the date of delivery of the equipment, the installation and training are free of charge, after exceeding this deadline, these services will be charged,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

