



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WSTECH MASTER PRO BOBATH Rehabilitation table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland, (2018),
- Upholstery with medical approval in green color has normal signs of use (possibility to replace the upholstery with a new one from our color chart for an additional fee),
- BOBATH is characterized by a stable powder-coated metal structure. The table is adjusted by means of a hydraulic actuator by pressing the lever located at the actuator several times with the foot,
- This model uses a bearing-based height adjustment mechanism. Thanks to the wide surface of the table, it is possible to work comfortably with children using the Vojta or Bobath method,
- Height adjustment by hydraulic actuator in the range of 48 - 100 cm,
- Overall dimensions: 110 x 200 x 48 / 100 cm,
- Maximum load capacity: 160 kg,
- Weight: 122 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

