



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DERUNGS HX35P FX Light (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 40 VA,
- Bulb: 6V / 35 W,
- Extremely lightweight design of the lamp allows for precise adjustment of the desired working position vertically and horizontally. Equipped with state-of-the-art LED light sources. Instant ignition and light stabilization thanks to the use of an electronic inverter,
- Elegant design and high quality and durability of the materials used. Thanks to its lightweight design, it is suitable for use both for examination in medical offices and as a friendly individual light source in patient rooms. The lamp mounts to a rail, wall, table or on a mobile tripod,
- FX version with rigid arm,
- Technical specifications:
 - Color rendering index: 100 Ra,
 - Color temperature: 3000 K,
 - LED lamp,
- Included:
 - Power cord,
- Dimensions: 10 x 15 x 105 cm,
- Weight: 1.3 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

