



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET ALPHASTAR 1132 Operating O/R table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: through 2 new built-in 12 V / 17 Ah batteries,
- Entire 5-segment table top transparent to X-rays. The table has REVERS function (inverted patient position),
- Frequency: 50/60 Hz,
- Power: 400 VA,
- Mobile universal operating table with electrohydraulic drive. Emergency control panel on the table column,
- New upholstery (memory visco) with medical approval in black,
- Adjustable in multiple planes by wired remote control or via panel on column:
 - ON/OFF,
 - LOCK/UNLOCK - lock and unlock the table,
 - NORMAL/REV - orientation of the patient on the operating table,
 - TREND/REV - change of the patient's lying direction in relation to the table column: Trendelenburg and anti-Trendelenburg positions within $\pm 30^\circ$,
 - The table is equipped with an auxiliary drive controlled by a hand-held remote control possible forward and backward travel,
 - UP/DOWN - operating surface height (up/down),
 - TILT L/TILT R - lateral tilt within 20° ,
 - BACK UP/BACK DN - back segment position adjustment,
 - LEG UP/LEG DN - adjustment of leg segment: up/down,
 - FLEX/REFLEX - bending the top in half in both directions: roof position / V position (Beach chair),
 - LEVEL - resetting: return to horizontal recliner position,
- Manually adjustable:
 - Headrest,
- Base housing made of durable plastic,
- The table features X - RAY translucency,
- Equipped with rails for attaching additional accessories,
- Table adapted for C-arm,
- Mobile, on double casters
- Manually removable headrest and foot segment,
- Includes:
 - Power cord,
 - Wired remote control,
- Possibility to retrofit the table at an additional cost with other accessories from our offer, examples shown in the last photos,
- Dimensions:
 - Headrest: 50 x 26 x 7 cm,
 - Backrest: 55 x 54 x 7 cm,
 - Seat: 55 x 42 x 7 cm,
 - Main leg rest segment: 55 x 52 x 7 cm,
 - Lower leg rest segment: 54 x 26 x 7 cm,
- Maximum load capacity: 450 kg,
- Weight: 325 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

