



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET ALPHASTAR 1132 Operating O/R table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real product photos,
- German production,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in batteries 12 V / 17 Ah,
- The entire 5-segment tabletop is radiolucent. The table features a REVERS function (inverted patient position),
- Frequency: 50/60 Hz,
- Power: 400 VA,
- Mobile universal operating table with electrohydraulic drive. Emergency control panel on the table column,
- New upholstery (memory visco) with medical certificate in black,
- Adjustable in many planes using a wired remote control or via the panel on the column:
 - ON/OFF,
 - LOCK/UNLOCK - locking and unlocking the table,
 - NORMAL/REV - patient orientation on the operating table,
 - TREND/REV - changing the patient's lying direction relative to the table column: Trendelenburg and anti-Trendelenburg position within the range of $\pm 30^\circ$,
 - The table is equipped with an auxiliary drive controlled by a handheld remote control, allowing forward and reverse travel,
 - UP/DOWN - height of the operating surface from 68.5 to 107 cm (up/down),
 - TILT L/TILT R - lateral tilt in the range of 20° ,
 - BACK UP/BACK DN - adjustment of the back section position,
 - LEG UP/LEG DN - adjusting the leg segment: up/down,
 - FLEX/REFLEX - bending the tabletop in half in both directions: roof position / V position (Beach chair),
 - LEVEL - reset: return to the horizontal position of the bed,
- Manually adjustable:
 - Headrest,
- The base housing is made of durable plastic,
- The table has X-RAY transparency,
- Equipped with rails for mounting additional equipment,
- Table adapted to C-arm,
- Mobile, on double wheels
- Manually removable headrest and foot section,
- Included in the set:
 - Power cord,
 - Wired remote control,
- The table can be equipped with other accessories from our offer for an additional fee, examples are shown in the last photos,
- Dimensions:
 - Headrest: 50 x 26 x 7 cm,
 - Backrest: 55 x 54 x 7 cm,
 - Seat: 55 x 42 x 7 cm,
 - Main leg rest segment: 55 x 52 x 7 cm,
 - Lower leg rest segment: 54 x 26 x 7 cm,
- Maximum load capacity: 450 kg,
- Weight: 325 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months to the international market,
 - Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

