



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ASPEL ASCARD RED3 ECG unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland, (2012),
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in Li-Ion 8.4 V / 2600 mAh battery,
- Frequency: 50/60 Hz,
- Power: 60 VA,
- AsCARD Red3 electrocardiograph is an innovative technical solution to ensure the highest quality and accuracy of examinations. Thanks to the applied touch panel, intuitive menu and clear function keypad, operation of the device is extremely easy and pleasant. The small compact size, low weight and battery power, further make the device mobile and convenient to use,
- Technical specifications:
 - Automatic analysis and interpretation,
 - Display of waveforms from 12 ECG leads,
 - Printing in 3, 6 or 12 channel mode,
 - Keypad with function keys,
 - On-screen menu for easy operation via touch panel,
 - Database of 10 recent examinations,
 - Perform up to 130 automatic examinations on battery operation,
 - Pacemaker detection,
 - Continuous measurement of heart rate (HR) and its presentation on the display,
 - Audible indication of heart rate,
 - Apparatus suitable for direct open-heart operation,
 - Automatic analysis and interpretation,
 - Network noise filter: 50Hz/60Hz,
 - Muscle interference filter; selectable filters: 25Hz, 35Hz, 45Hz,
 - Isoline filter; selectable filters: 0.15Hz, 0.45Hz, 0.75Hz, 1.5Hz,
 - Electrode detachment detection independent for each channel,
 - Print accuracy setting (line thickness of printed curves),
 - USB interface for communication with PenDrives, allowing directly from the camera to record examinations on a PenDrive in CardioTEKA standard and its subsequent automatic analysis and interpretation in CardioTEKA software,
 - Viewing on the display of the examination descriptions stored in the memory, with the possibility of editing the data commenting on the examination,
- Selectable menu in language: PL, EN,
- Keyboard on the control panel in language: English,
- Included:
 - User manual in Polish (PDF),
 - EKK clip limb electrodes - 4 pcs,
 - EPP suction cup electrodes - 6 pcs,
 - ECG 10-lead cable,
 - Power cord,
- Dimensions: 28 x 47 x 12 cm,
- Weight: 3 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

