



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## OLYMPUS OTV-S7H 1D L08E Camera head (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- OLYMPUS OTV-S7 camera heads are high quality medical equipment, which is used daily by doctors around the world, helping thousands of patients. The presented model OTV-S7H-1D-L08E was created by OLYMPUS using modern technologies,
- The camera head of OTV-S7 uses 1080i resolution to provide a sharp and accurate image, which is essential in surgical procedures,
- Endoscope input:  $\varnothing 45.2 \times 12.1$  mm,
- Disinfection / sterilization:
  - Ethyl oxide sterilization capability,
  - Can be disinfected in disinfectant solution,
- Dimensions:
  - Camera head:  $\varnothing 14 \times 45$  mm,
  - Cable length: 4 m,
- Weight: 0.4 kg,
- Has an up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

