



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS OTV-S6H-NA Camera head (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- OLYMPUS OTV-S6H-NA is a high quality medical equipment, which is used daily by doctors around the world, helping thousands of patients. The presented model OTV-S6H-NA was created by OLYMPUS using modern technologies,
- Built-in 2 remote control buttons,
- Retrofitted OLYMPUS ART-T12QA adapter,
- Disinfection / sterilization:
 - Ethyl oxide sterilization capability,
 - Possibility of sterilization by steam (autoclave),
 - Possibility to disinfect the adapter with disinfectant solution,
- Included:
 - OLYMPUS container for sterilization,
- Dimensions:
 - Cable length: 4 m,
- Weight: 0.3 KG,
- Has up-to-date inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

