



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Medical kidney dish (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- Real photos of the product,
- Metal kidney dish is a reusable kidney-shaped medical bowl. It belongs to the basic equipment of medical offices, emergency team or hospitals,
- The kidney-shaped tray is made of extremely resistant non-magnetic stainless steel with parameters that meet high quality standards,
- Specifications:
  - High quality workmanship,
  - Non-magnetic stainless steel,
  - Smooth edge ends,
- Dimensions: 24 x 12 x 4.5 cm,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

