



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Rigid gynecological scraper (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Gynecological scraper has blunt-rigid blades, which minimizes the risk of damage and irritation during procedures, while ensuring precise and effective operation,
- Ergonomic handle guarantees a firm grip and control of the tool, resulting in safe and effective procedures,
- Realistic product photos,
- Manufactured from stainless steel,
- Satin finish,
- Fully autoclavable scraper,
- Neck length: 19 cm,
- Has an up-to-date inspection and is ready for operation,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

