



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MOBICLINIC QDC-303 Anti-decubitus mattress (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Spain, (2020),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- The anti-decubitus bubble mattress MOBICLINIC QDC-303 is an effective prevention and supports the treatment of decubitus ulcers up to and including the first degree. Thanks to the alternating inflation of the chambers, the so-called "continuous motion effect", which provides changes in pressure on different areas of the body, as well as provides continuous muscle massage and improves blood flow,
- This makes the mattress an excellent tool for the prevention of bedsores. It provides pressure relief to the sensitive areas most prone to pressure sores, such as the heels, elbows, buttocks, shoulder blades, occiput,
- Features:
  - Infinite manual pressure adjustment of the mattress, individually selected according to the patient's weight,
  - Quiet energy-saving pump - important information when the pump must work all the time, including night,
  - The pump is equipped with hooks that allow vertical fixation,
  - Thanks to the design based on independent air chambers, even distribution of the patient's weight, proper ventilation and, thanks to the excellent insulating properties of the air, temperature comfort is achieved,
  - Work cycle: 6 minutes,
  - Pads for attachment to regular bed mattresses,
- Technical specifications:
  - Pressure range: 40 - 100 mmHg,
  - Air flow rate: 7-8 l/min,
- Includes:
  - User's manual in language (PDF),
  - Power cord,
- Dimensions:
  - Pumps: 13.5 x 26.5 x 10 cm,
  - Mattress: 90 x 200 x 7 cm,
- Maximum load capacity of the mattress: 135 kg,
- Weight: 7 kg,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

