



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



PENTAX EC-3885LK Videocolonoscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Colonoscope EC-3885LK provide excellent high-resolution imaging in a full line of economical probes with ergonomic design and established quality,
- Includes high-resolution imaging sensors with bright illumination to ensure effective image clarity and clear details,
- Technical specifications:
 - Working length of the viewing section: 1700 mm,
 - Working part diameter: 12.8 mm,
 - Diameter of operating channel: 4.2 mm,
 - Viewing direction: forward,
- Maximum tip bending angle:
 - Top: 180°,
 - Bottom: 180°,
 - Right: 160°,
 - Left: 160°,
 - Viewing angle: 140°,
- Compatible with:
 - EPK-700,
 - EPK-1000,
 - EPK-i,
 - EPK-i5000,
 - EPK-i7000,
- Includes:
 - Pentax case,
 - Cover,
- Dimensions: 16 x 59 x 43 cm,
- Weight: 4 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

