



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER etCO2 module (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Frequency: 50/60 Hz,
- The Drager etCO2 is an advanced module designed to accurately monitor carbon dioxide (CO2) in the respiratory system of a patient undergoing mechanical ventilation. Integrated with Dräger masks, this module provides continuous and timely patient monitoring, enabling accurate assessment of lung function and early detection of respiratory complications,
- Specifications:
 - Continuous CO2 monitoring:
 - The CO2 module monitors carbon dioxide exhalation by monitoring the precise assessment of alveolar ventilation and lung function,
 - Precision CO2 sensor:
 - Equipped with a high-precision CO2 sensor, the module provides precise and reliable measurement of CO2 concentration, allowing early detection of respiratory abnormalities,
 - Configurable alarms:
 - The module is equipped with configurable alarms that alert physicians to abnormal CO2 concentrations, allowing quick intervention,
 - Real-time display:
 - CO2 levels are displayed at regular intervals on the breathing apparatus, providing continuous monitoring and immediate assessment of respiratory function,
 - Data storage:
 - The module records and stores CO2 data for further analysis and accurate documentation of patient ventilation,
- Advantages:
 - Continuous and precise monitoring of CO2 levels for accurate assessment of lung function,
 - Early detection of respiratory complications and rapid intervention in case of injury,
 - Increased patient safety with customizable alarms and immediate monitoring,
- Includes:
 - ZOLL CAPNOSTAT sensor,
 - Monitoring cable,
- Dimensions: 25 x 25 x 5 cm,
- Weight: 0.5 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

