



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



GRAND SEIKO GR-2100 Autorefractometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 80 VA,
- LCD color display with a screen diagonal of 5.6",
- Has the following diagnostic capabilities:
 - Measurement modes: Auto / Manual / Quick,
 - IOL Option,
 - Patient VD: 0, 10, 12, 13.5, 15 mm,
 - Sphere (SPH): ± 25 D,
 - Cylinders (CYL): ± 10 D,
 - Form of cylinders: +/-mix,
 - Axis (Axis): from 10 to 180°,
 - Minimum pupil diameter: 2.3 mm,
- It has a built-in thermal printer with automatic paper cutting, paper width: 57 mm,
- Power saving mode: Off after 3, 5 or 10 minutes,
- Manual joystick canopy adjustment of 38 mm,
- Manual chin adjustment within 60 mm,
- Output port: RS232C, Video-out,
- Dimensions: 27 x 42 x 42.5 cm,
- Weight: 15 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

