



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DATEX-OHMEDA TEC 7 Sevoflurane vaporizer (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- DATEX-OHMEDA TEC 7 vaporizers allow efficient and economical administration of inhalant,
- Provides a constant initial vapor concentration over the entire range of flow rates used in clinical practice,
- Easy-to-turn dial and fine adjustment capability helps to precisely control the concentration of the administered anesthetic,
- The vaporizer knob for sevoflurane allows step changes in 0.2% increments for concentrations up to 1% and 1% increments for concentrations from 1% to 8%,
- To reduce the frequency of filling, the TEC 7 vaporizer holds 225 ml of anesthetic,
- The TEC 7 vaporizer can be installed in Datex-Ohmeda anesthesia machines equipped with the Selectatec attachment manifold. No tools are required for this operation, and the anesthesia machine does not need to be taken out of service,
- TEC 7 vaporizers can be installed on a single mounting manifold with Datex-Ohmeda TEC 4, TEC 5, TEC 6 and TEC 6 Plus vaporizers. The listed vaporizers are compatible with the co-blocking mechanism that the TEC 7 vaporizer is also equipped with,
- Ability to read the level of liquid anesthetic thanks to an indicator on the front of the evaporator,
- Selectatec (S - 2000) attachment is used in the given vaporizer,
- Fill type: Quick Fill,
- Included:
  - English instruction manual (PDF),
- Dimensions: 22 x 11 x 25 cm,
- Weight: 8 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

