

PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FRESENIUS KABI INJECTOMAT MC Syringe infusion pump (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- French production,
- Power supply: 230 V,
- Emergency power supply: new built-in battery 6V/1.8Ah,
- Frequency: 50/60 Hz,
- Power: 15 W,
- **FRESENIUS KABI INJECTOMAT MC** single-syringe pump, electronically controlled, enabling cooperation with the central power supply and data management system,
- The pump is calibrated to work with syringes of various types and volumes and from various manufacturers,
- **Characteristics:**
 - Versatile syringe pump with dose calculation,
 - To conduct pharmacological treatment,
 - Advanced features include dose calculation in mass units, programmable bolus, initial dose,
 - Easy to use,
 - Intuitive operation,
 - Clear, interactive menu,
 - Durable and proven construction,
 - Library of 19 drugs,
- **Technical specifications:**
 - Range, dosing rate in the range of 0.1-1200 ml/h for a 50/60 ml syringe.,
 - BOLUS dose rate in the range of 0.1-1200 ml/h for a 50/60 ml syringe,
 - Pressure modes: 2 modes available: variable or 3 pre-set levels - from 100 to 900 mmHg, in 50 mmHg increments,
 - DPS Dynamic Pressure System (DPS) alerts you to pressure changes. This allows you to predict the risk of blockage or possible leakage in the infusion line.
 - Anti-bolus system: Significantly reduces the bolus after releasing the occlusion (max 0.2 ml),
 - Possibility to change the speed without stopping the infusion and to administer a bolus dose at any time during the infusion at any speed,
 - Blue, monochrome graphic LCD display, 70 mm x 35 mm,
 - "Night mode" with the option of switching between night and day mode manually and automatically,
 - Possibility of programming infusion parameters in various units, taking into account body weight in relation to time: mg/ng/kg/g/mmol/mU/U/kcal/mEq,
 - Possibility to connect pumps into modules without using a docking station,
 - Drugs contained in the drug library linked to infusion parameters, ability to display the drug name and/or selected infusion parameters,
- **The pump supports syringes:**
 - ASTRA ZENECA 50 ml.,
 - BD PERFUSOR 50 ml.,
 - BD PLASTIPAK 50/30/20/10/5 ml.,
 - BD PLASTIPAK 50 WWD 50/20 ml.,
 - BRAUN OMNIFIX 50 ml.,
 - BRAUN PERFUSOR 50/20 ml.,
 - CODAN 50/10 ml.,
 - DIDACTIC LINE 50 ml.,
 - DIDACTIC PERFUSION 50 ml.,
 - DISPOMED INFUSECT 50 ml.,
 - DISPOMED PERFUSECT 50 ml.,
 - FRESENIUS INJECTOMAT 50/10 ml.,
 - FRESENIUS P-SPRITILE 50 ml.,
 - ICOGAMMA PLUS 50 ml.,
 - SANPOL 50 ml.,
 - KABIFILL 50 ml.,
 - MAROOMED 50 ml.,
 - MOWOSECT 50/30/20/10/5 ml.,
 - PENTA 60/30/20/10 ml.,
 - PIC INDOLOR 50 ml.,
 - POLFA 50/20/5 ml.,
 - TERUMO 50/30/20/10/5 ml.,
 - TUTOJECT TYPET 50 ml.,
 - BRAUN OMNIFIX 30/20 ml.,
 - FRESENIUS MED.CARE 30/20 ml.,
- **Possibility to select menu in: PL, FR, EN, DE, NR, SE, DU, IT, ES, GR, RU, HU, CS, SK, TR, RU, PT, BR,**
- **Included in the set:**
 - User manual in German (PDF),
 - Power cord,
- Dimensions: 13.5 x 34.5 x 17 cm,
- Weight: 2.1 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for a period of 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

