



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZEISS IOL MASTER Optical biometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 90 VA,
- Software version: 5.0.4.,
- Technical specifications:
 - Axial length:
 - Area: 14 - 40 mm,
 - Display scaling: 0.01 mm,
 - Keratometer:
 - Area: 5 - 10 mm,
 - Display scaling: 0.01mm,
 - Anterior chamber depth:
 - Area: 1.5 - 6.5mm,
 - Display scaling: 0.01mm,
 - White off white:
 - Area: 8 - 16mm,
 - Display scaling: 0.1 mm,
 - Adjustment of chin rest in range: 20 - 26 cm,
- Selectable menu in the following languages: DE, EN, FR, ES, IT, SE,
- Included:
 - HP DESKJET 5150 printer,
 - Wired keyboard and mouse,
 - IOL Master Artificial Eye,
 - CARL ZEISS mobile table with electric height adjustment from: 64 - 92 cm,
 - User manual in English (PDF),
 - Power cords,
- Dimensions:
 - Tabletop: 90 x 42 cm,
 - Biometer: 32 x 40 x 58 cm,
 - Overall: 90 x 55 x 123 cm,
- Total weight: 55 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

