



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **MGE SAM 14 Surgical suction unit (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the UK, (2019),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 90 W,
- The SAM 14 surgical suction machine was designed for quiet and efficient operation through an oil-free diaphragm pump,
- The SAM 14 Respirator is a top-quality surgical suction machine designed for high customization and a high-flow, high-vacuum device ideal for use in general wards, emergency rooms, smaller operating rooms, intensive care, and as a backup device in hospitals and dental clinics,
- Technical specifications:
  - Airflow: 30l / min,
  - Vacuum: 0 to 0.85 bar / 0 to 650 mmHg,
- Included:
  - Power cord,
  - 2 pcs. - 2 l jar. SERRS,
  - 2 pcs. disposable jar cartridge,
  - 2 pcs. - biological filter,
- Dimensions: 44 x 24 x 74 cm,
- Weight: 10 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



