



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS S4-2 Sector head (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the PRC,
- Technical specifications:
 - Frequency range: 2 - MHz,
 - Field of view: 90°,
 - Number of elements: 80,
 - Supported modes: 2D, CW, pulse wave, color Doppler, tissue Doppler, XRES, AutoSCAN, iSCAN and harmonic imaging,
 - Designed for use in adult, abdominal, pediatric and TCD echocardiography,
 - Probe type:
 - Sectoral,
 - Broadband,
- Dimensions: 8 x 3.5 x 2 cm,
- Weight: 0.5 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

