



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZOLL M SERIES Biphasic defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 10 - 29 V,
- Emergency power supply: via 1 new built-in 10 V / 2.5 Ah battery,
- Frequency: 50/60 Hz,
- Power: 130 W,
- The Zoll M Series is a lightweight and portable device that combines defibrillation, ECG, patient monitoring functions with the ability to print and record data,
- The defibrillator is designed for use during any type of resuscitation, and its light weight, small size and compact design make it particularly suitable for use during transport,
- It has:
 - 5.5" color LCD screen,
 - Menu in language: English,
 - Keypad on the control panel in: English,
 - Built-in module:
 - ECG - 10-lead ECG cable,
 - SPO2,
 - RS-232/EKG X 1000 connector,
 - Printer for 90 mm wide paper,
 - Port for external PCMCIA memory cards,
 - Portable thanks to convenient handle and wiring bag,
- Includes:
 - User manual in Polish (PDF),
 - Power cord,
 - Dual-function (adult-child) hard spoons,
 - SPO2 cable,
 - 10-lead ECG cable,
- Dimensions: 36 x 47 x 22 cm,
- Weight: 10 kg,
- Cat D136,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

