



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **DRAGER FABIUS CE Anesthesia machine (Respirator/Ventilator) (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 3.2 Ah batteries,
- Frequency: 50/60 Hz,
- Software: 02.80,
- Has two ports for vaporizers with common Dräger DW-2000 connector systems, no vaporizers included - can be retrofitted at additional cost,
- Respirator driven by an electric motor (electronically controlled),
- Airway monitor with control unit is a device for monitoring and measuring ventilation parameters during anesthesia,
- The monitor is equipped with a control unit to set and change all ventilation parameters listed below,
- Supplied with gases:
  - 2 x O<sub>2</sub>,
  - 2 x N<sub>2</sub>O,
  - 2 x Air,
- Ventilator measurement functions:
  - Maximum pressure (P<sub>max</sub>) 10 to 70 hPa,
  - Tidal volume (VT) 50 to 1400 mL,
  - Respiratory rate (f<sub>IPPV</sub>) 6 to 60 L/min,
  - Inspiration/Expiration Ratio, (Ti : Te) 1 :3 to 2: 1,
  - Inspiratory pause, (Tip : Ti) 5 to 50%,
  - Positive end expiratory pressure (PEEP) 0 to 15 hPa,
- Infinitely variable manual adjustment of gas flow - O<sub>2</sub>, AIR, N<sub>2</sub>O, 0.0 to 12L/min, (pressure gauges),
- Safety functions:
  - P, V, FiO<sub>2</sub> monitoring,
  - Automatic O<sub>2</sub> shortage alarm and N<sub>2</sub>O cutoff,
  - It has a built-in S-ORC (Sensitive Oxygen Ratio Controller), a Sensitive Oxygen Ratio Controller (a control device to ensure a minimum O<sub>2</sub> concentration of 21 Vol. %),
- Compact patient breathing system with oxygen sensor, located on a movable arm,
- Built-in electronic display with O<sub>2</sub>, AIR, N<sub>2</sub>O gas measurement function, and control panel and ventilator function display,
- Clear and easy to use and interface,
- Equipment:
  - 3 drawers with dimensions: 32 x 30 x 11.5 cm,
  - Slats for attaching additional accessories,
  - 2 holders for gas cylinders,
- Mobile machine on wheels, including 2 with brakes,
- Selectable menus in ENG, GER, FR, ESP, IT,
- Keypad on control panel in language: English,
- Apparatus equipped with:
  - HEWLETT PACKARD VIRIDIA 24C Patient monitor / cardiac monitor:
    - Reconditioned (used),
    - Technical condition: very good,
    - Visual condition: very good,
    - Made in the USA,
    - Power supply: 230 V,
    - Frequency: 50/60 Hz,
    - Technical specifications:
      - Liquid crystal monitor with screen diagonal: 11",
      - Foiled and easy to read control panel,
      - Audible and visual alarms,
      - Equipped with a docking station following with modules:
        - SPO<sub>2</sub>,
        - NiBP,
        - ECG,
  - HEWLETT PACKARD gas measurement module:
    - Automatic alarm when anesthetic mixture is detected,
    - Inspiratory and exhaled values for all gases,
    - Fast warm-up time,
    - Low sample flow rate,
    - Individual alarm thresholds for each parameter,
    - Automatic timer calibration,
    - Respiratory rate calculated based on CO<sub>2</sub> waveform,
- Included:
  - Gas supply hoses with tips of customer's choice (AGA or DIN),
  - User manual of the DRAGER FABIUS CE device in Polish (PDF),
  - Power cord,
  - New SPO<sub>2</sub> clip sensor,
  - New 3-lead ECG cable,
  - NiBP cuff with air tubing,
- Dimensions: 70 x 102 x 154 cm,
- Weight: 113 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - For the domestic market (Poland), if the device is installed within 30 days from the date of delivery of the equipment, installation and training are free of charge, after exceeding this period, these services will be charged,
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty to 12-months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

