



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **DRAGER D-VAPOR 3000 (Desofluran/Desoflurane) inhalation agent vaporizer (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production, (2013),
- Power supply: 230 V,
- Memory backup power supply: via 1 new built-in 6V battery,
- Frequency: 50/60 Hz,
- Software: 3.00,
- **The D-Vapor 3000** is an anesthetic vaporizer designed for use with Drager anesthesia machines,
- It is used to vaporize desflurane, a strong inhalation anesthetic,
- **The D-Vapor 3000** is a precise and reliable vaporizer that offers many features that improve safety and ease of use,
- **Characteristics:**
  - The volume of the vaporizer is 300 ml, so the entire contents of a standard 240 ml anesthetic bottle can be poured into the vaporizer, even if there is anesthetic left inside from the previous filling,
  - Unique transport feature - vaporizers can be carried in any position because the transport lock prevents the anesthetic from spilling, so it can be used immediately after transport,
  - All filling systems are designed specifically for specific inhalation anesthetics, eliminating the risk of accidentally using the wrong anesthetic.
- **Technical specifications:**
  - Flow range: From 0.2 to 15.00 (L/min),
  - Possibility to adjust the desoflurane concentration in the range of 2 - 18%,
  - Infusion type: Drager Fill,
  - This evaporator uses an AUTO-EXCLUSION mounting,
- **Included in the set:**
  - Power cord,
- Dimensions: 11 x 26 x 26 cm,
- Weight: 6.4 kg,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

