



DRAGER VAPOR 19.3 Sevoflurane vaporizer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Vapor 19.3 is an anesthetic vaporizer, allows efficient and economical administration of inhalant,
- Technical specifications:
 - Vaporizer volume: max 200 ml,
 - Operating temperature: 15° C - 35° C,
 - Infinitely adjustable anesthetic concentration,
 - Unique transport function - vaporizers can be carried in any position. The transport safety device prevents the anesthetic from spilling, so it can be used immediately after transporting,
 - All filling systems are designed specifically for specific inhalation anesthetics. This means eliminating the risk of mistakenly using the wrong anesthetic,
 - The particular vaporizer uses a Dräger mount (DW-2000),
 - Fill type: Funnel Fill,
- Dimensions: 10 x 17 x 24 cm,
- Weight: 8 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

