



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybierański
+48 882 762 006

kw@praiston.pl



ASCOR AP12 Infusion (syringe) pump (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: via new built-in 9.6 V / 2000 mA battery,
- Frequency: 50/60 Hz,
- Power: 10 VA,
- The AP12 single syringe pump is designed for use in anesthesia, intensive care, cardiology, neonatology, oncology,
- The AP12 pump has a very user-friendly dose programming system with the ability to set: speed, speed and volume, speed and time, or volume and time,
- Features:
 - Ease of use and reliability,
 - Ability to view and change parameters during infusion,
 - Working with 10-60 ml syringes (about 40 types),
 - BOLUS safe shock dose delivery function,
 - Extensive alarm system,
 - Built-in drug library,
 - Long-term battery operation possible,
 - Built-in test system,
 - Can be mounted on a tripod, bed or other rail,
- Technical specifications:
 - Dispensing speed:
 - Adjustable in 0.1 ml/h increments from 0.1 ÷ 99.9 ml/h and in 1 ml/h increments from 100-500 ml/h:
 - 0.1 ÷ 500 ml/h for 50 ml syringe,
 - 0.1 ÷ 300 ml/h for 30 ml syringe,
 - 0.1 ÷ 250 ml/h for 20 ml syringe,
 - 0.1 ÷ 150 ml/h for a 10 ml syringe,
 - Maximum rate of shock dose (bolus) dosing:
 - Adjustable in 1 ml increments:
 - 2000 ml/h for 50 ml syringe,
 - 1200 ml/h for 30 ml syringe,
 - 1000 ml/h for 20 ml syringe,
 - 600 ml/h for a 10 ml syringe,
 - Shock dose dispensing volume:
 - 0.1 ml to the volume of the selected syringe, programmable in 0.1 ml increments,
 - KVO maintenance dose:
 - 0 ÷ 5 ml/h, programmable in 0.1ml/h increments,
 - Dosage accuracy:
 - ± 2% (according to the requirements of EN 60601-2-24),
 - Infusion volume:
 - 0.1 ÷ 1000ml, programmable in 0.1ml and 1ml increments from 100 to 1000ml,
 - Infusion time:
 - Max 99 hr 59 min,
 - Occlusion pressure:
 - Minimum: 0.04 MPa ±0.01 (300 mm Hg ± 75),
 - Low: 0.06 MPa ±0.02 (450 mm Hg ±150),
 - Medium: 0.09 MPa ±0.02 (675 mm Hg ±150),
 - High: 0.12 MPa ±0.03 (900 mm Hg ±225),
 - Syringe types:
 - BRAUN OMNIFIX 10 ml,
 - BOGMARK 10 ml,
 - CODAN 10 ml,
 - ERG 10 ml,
 - MONOJECT 10 ml,
 - TERUMO 10 ml,
 - BRAUN OMNIFIX 20 ml,
 - BRAUN PERFUS 20 ml,
 - BD PLASTIPAK 20 ml,
 - BOGMARK 20 ml ..,
 - CODAN 20 ml ..,
 - ERG 20 ml ..,
 - MONOJECT 20 ml,
 - POLFA LUBLIN 20 ml ,
 - SCHAN CHUAN Syringe Set 20 ml./cc,
 - BRAUN OMNIFIX 30 ml,
 - BD PLASTIPAK 30 ml,
 - ONCE 30 ml,
 - MONOJECT 30 ml,
 - TERUMO 30 ml,
 - BD PERFUSION 50 ml,
 - BD PLASTIPAK 50 ml,
 - BRAUN ADALAT 50 ml,
 - BRAUN OMNIFIX 50 ml,
 - BRAUN PERFUSOR 50 ml,
 - PROINJEKT Spritze B Braun 50 ml,
 - CODAN/ONCE 50 ml,
 - DISPOMED 50 ml,
 - DISPOMED type P 50 ml,
 - DISPOMED 50 ml,
 - INJECTOMAT 50 ml,
 - IVAC 50 ml,
 - JANPOL 50 ml,
 - MONOJECT 50 ml,
 - ONCE PERFUSION 50 ml,
 - POLFA BOLESŁ. 50 ml,
 - POLFA BOLESŁ. 50 ml,
 - POLFA LUBLIN 50 ml,
 - SCHAN CHUAN Syringe Set 50 ml./cc
 - TERUMO 50 ml,
 - PENTA (PF) 60 ml,
- Includes:
 - English user manual (PDF),
 - Power cord,
- Dimensions: 18 x 32 x 14 cm,
- Weight: 3.1 KG,
- Cat 02,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

