



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WALDMAN UV181 BL-UV200 L Therapy light (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- This model is designed for the examination of small sections of the limbs. The basic UV181 model is designed to treat the hands and lower extremities. It also has a built-in ventilation device to minimize surface temperature and, like the UV200AL, is used to treat the hands and arms,
- By partially irradiating the body with the UV181 and UV200, conditions of the hands (such as atopic palm eczema) and soles of the feet can be effectively treated,
- UV181 with adjustable horizontal tilt in the range: 0 / 90°,
- With the built-in LCD timer, doses can be entered in joules per cm² or minutes/seconds,
- Included:
 - Power cord,
- Dimensions: 64 x 70 x 125 cm,
- Weight: 83 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

