



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE ICC-200 Surgical diathermy (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 200 W,
- Technical specifications:
 - Infinitely adjustable power,
 - Digital displays of set power,
 - Foiled, touch-sensitive control panel,
- Features the following functions:
 - ENDO-CUT,
 - Autocut - max. 200 W,
 - Autocoagulation - max. 120 W,
 - Available coagulation modes: SOFT, FORCED, BIPOLAR,
- Equipment:
 - Foot switch,
 - Brand new bipolar cord, BOWA forceps (bipolar tweezers can be retrofitted at additional cost),
 - Brand new monopolar electrode lead with BOWA grip - two buttons, 4 mm shaft (no electrodes included, can be retrofitted at additional cost),
 - Brand new cable with BOWA clip for disposable neutral electrodes (no electrodes included, can be retrofitted at additional cost),
- Keypad on control panel in language: English,
- Included:
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 28 x 38 x 15 cm,
- Weight: 8.4 KG,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

