



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDIFA MAT 5000 Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2012),
- MEDIFA MAT 5000 is designed to perform procedures and surgeries in the field of general surgery. All metal parts exposed to external factors are made of acid-resistant steel. The change of the height of the table top is possible thanks to a hydraulic actuator and foot pump and is carried out with a foot switch. The base of the table is mobile, with central wheel lock and table positioning function activated by a foot switch. Changing the angle of the backrest, footrest, headrest, as well as the Trendelenburg and anti-Trendelenburg positions is possible with the use of interlocking gas springs, actuated by levers located in convenient places for the operator. Lateral tilting is accomplished with a crank located on the side, under the seat segment. The tabletop segments are permeable to x-rays and allow imaging and C-arm monitoring. The limitation of the monitoring area is the width of the table column. The table mattresses are seamless, ergonomically shaped and transparent to X-rays,
- Upholstery with medical approval in black color,
- 4 - segmented,
- Manually adjustable in many planes:
  - Height in range: 74 - 99 cm,
  - Headrest angle in range: from -45° to +45°,
  - Backrest angle in the range of -40° to +70°,
  - Lateral tilt:  $\pm 20^\circ$ ,
  - Trendelenburg position: 30°,
  - Anti-Trendelenburg position: 30°,
  - Leg section in range: from -90° to +25°,
- Features:
  - X-Ray translucency,
  - Rails for attaching additional accessories,
  - Stainless steel base,
  - Mobile with brake,
- Dimensions:
  - Base: 60 x 100 cm,
  - Bed: 54 x 205 cm,
  - Headrest: 50 x 30 cm,
  - Backrest: 54 x 62 cm,
  - Seat: 54 x 43 cm,
  - Leg rest: 54 x 63 cm,
- Maximum load capacity: 150 KG,
- Weight: 130 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

