



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDTRONIC BIPHASIC LIFEPAK 12 Defibrillator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: via 2 built-in 12 V / 1.2 Ah batteries (one new included, can be retrofitted with a second new one at additional cost),
- Frequency: 50/60 Hz,
- Power: 180 W,
- Selectable menu in the following languages: EN, FR, DE, ES, SE, IT, FI, DU, PL, PT, DE, NR,
- Keypad on control panel in language: English,
- Built-in printer for paper width: 11 cm,
- Equipped with a carrying handle,
- Features:
  - SYNC - Synchronization - Cardioversion,
  - ADVISORY - INFO (advisory) - Activates the continuous patient surveillance system (CPSS),
  - ANALYZE - ANALYZE - Activates the SAS (Shock Advisory System),
  - PACER - Stimulator - Activates the stimulation function,
- Equipped with the following modules with wiring and sensors:
  - ECG - 10-lead cable,
  - NiBP cuff,
  - SPO2 sensor,
  - TEST LOAD Tester,
  - QUICK COMBO cable,
- Defibrillator equipped with a dedicated bag with side pockets for transporting necessary instrumentation,
- Includes:
  - User manual in Polish (PDF),
- For an additional fee, the defibrillator can be retrofitted with:
  - An external charger from AD Elektronik shown in the last photo,
  - Additional battery,
  - Hard spoons,
- Overall dimensions: 60 x 20 x 33 cm,
- Weight: 10 kg,
- Kat 79,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

