



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER OXYLOG 2000 Ambulance ventilator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 12 V,
- Emergency power supply: via new built-in 1000 mAh battery,
- Frequency: 50/60 Hz,
- Technical specifications:
 - Built-in pressure gauge: -10 to 80 mbar,
 - Ventilation frequency (Freq): 0 - 40 1/min,
 - Tidal volume (Vt): 0,1 - 1,5 L,
 - Ratio of inspiratory and expiratory time (Ti:Te): 1:3 / 1:2 / 1:1,5 / 1:1 / 2:1,
 - Maximum airway pressure (Pmax): 20 - 80 mbar,
 - Pressure at the end of expiration (PEEP): 0 - 15 mbar,
 - Ventilation mode switch: SIMV/CPAP or IPPV/SIPPV,
 - Mixing switch: 100% O2 or 60% O2,
- Includes:
 - User manual in Polish (PDF),
 - AC adapter,
 - Oxygen supply hose with tips of the customer's choice (DIN or AGA),
 - Air flow valve,
 - O2 hose,
 - Can be retrofitted at additional cost:
 - Breathing valve (valve), Flow sensor, Angled connector,
 - Patient breathing tube set,
 - Oxygen cylinder 2 L,
 - Oxygen regulator,
- Dimensions: 21.5 x 21 x 17 cm,
- Weight: 7.1 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

