



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



DRAGER OXYLOG 2000 Ambulance ventilator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 12 / 30 V,
- Emergency power supply: through a new built-in 1000 mAh battery,
- Frequency: 50/60 Hz,
- Technical specifications:
 - Built-in pressure gauge: -10 to 80 mbar,
 - Ventilation frequency (Freq): 0 - 40 1/min,
 - Tidal volume (Vt): 0,1 - 1,5 L,
 - Ratio of inspiratory and expiratory time (Ti:Te): 1:3 / 1:2 / 1:1,5 / 1:1 / 2:1,
 - Maximum airway pressure (Pmax): 20 - 80 mbar,
 - Pressure at the end of expiration (PEEP): 0 - 15 mbar,
 - Ventilation mode switch: SIMV/CPAP or IPPV/SIPPV,
 - Mixing switch: 100% O2 or 60% O2,
 - Instrumentation housed in a dedicated carriage mount with handle,
- Includes:
 - User manual in Polish (PDF),
 - AC adapter with car adapter,
 - Oxygen supply hose with choice of tip (DIN or AGA),
 - Breathing valve, flow sensor, angle connector,
 - Breathing tube kit,
 - Air flow valve,
 - O2 hose,
- Can be retrofitted with 2 L. oxygen cylinder and oxygen regulator at additional cost,
- Overall dimensions: 50 x 27 x 53 cm,
- Weight: 11.5 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

