



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LID Julia LED light (New)**

- New,
- Made in France,
- Julia is the lamp of unique design with two powerful light beams (10° and 25°), which can be alternately switched with a two-position switch,
- The 10° beam provides strong, concentrated and precise illumination, ideal for local examination,
- The 25° beam provides diffuse and uniform illumination when examining a large area,
- Ergonomic handle provides fast and precise positioning under all circumstances,
- Specifications:
  - Light intensity: 740 Lumens,
  - Light color: 4500°K\* (neutral white),
  - Color rendering index: IRC 95/100,
  - Electricity consumption: 6,5 W,
  - Class 1 medical device,
  - LED lifetime: minimum 50,000 hours,
  - Electrical protection: insulation class II,
  - Standards: IEC 60 601 -1/ -2-41 IEC 60 598-2-25,
  - Power cord length: 3 meters,
  - Available colors : white,
- Dimensions:
  - Lamp head dimension: 100 x 110 mm,
  - Flexible arm: 650 mm, 1140 mm,
  - Base dimension: 69 x 69 x 115 mm,
  - Weight: 1.5 / 2.4 kg,
  - Inter-axis spacing (4 holes in lamp base): 54.5 mm,
- The standard length of the lamp is 650 mm, for an extra charge you can purchase a lamp with a length of 1140 mm
- Light field diameter :
  - Large beam 25 °:
    - 30 cm distance: 30,000 lux illumination field 15 cm,
    - Distance 50 cm: 11,000 lux field illumination 25 cm,
    - Distance 100 cm: 3,000 lux field illumination 35 cm,
  - Small beam 10 °:
    - Distance 30 cm: 110,000 lux field illumination 8 cm
    - Distance 50 cm: 45,000 lux field illumination 15 cm,
    - Distance 100 cm: 12,000 lux field illumination 25 cm,
- Additional mounting accessories are required for the lamp , of your choice (at an extra charge):
  - Telescopic stand base 60-100 cm white metal, heavy,
  - Telescopic stand base 60-100 cm white metal, heavy -version with brakes in 2 wheels,
  - Telescopic stand base 90-130 cm metal white, heavy,
  - Telescopic stand base 90-130 cm metal white, heavy -version with brakes in 2 wheels,
  - Plain clamp - black,
  - Clamp for medical rails with vertical lamp outlet,
  - Universal clamp,
  - Wall bracket,
  - ARM wall bracket white (2x40cm),
  - ARM wall mount white (1x40cm),
  - Tabletop mount ARM white (1x60cm),
  - Metal base white / gray,
  - Plate with metal pipe in the middle ø16,
  - Plate with metal tube on the edge ø16,
  - Plate on ARM holder,
  - Female socket with adaptor plate Ø 16 mm x 27,
- It has a current Technical Passport issued,
- Warranty: 3 years,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

