



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG Careclave 618 Autoclave (New)

- New,
- Made in Germany,
- With Careclave it is possible to clean and sterilize up to 8 dental handpieces, prophylactic handpieces or tips in as little as 22 minutes,
- Careclave cares for handpieces and turbines gently and quickly. In addition, the integrated ADDcare care station also makes rotor care directly with Careclave,
- Careclave uses patented DRYtelligence technology and a double jacket, making it offer everything expected of a modern autoclave. In addition, it enables sterilization of loads of up to 11 kg in up to 8 trays or 4 MELAstore containers in a short time. Based on the revolutionary chamber concept, washing, care and sterilization of instruments in the Carebox and simultaneous sterilization of other instruments in short trays is possible in one cycle,
- The ProControl system puts an end to the era of paper documentation. Now you can quickly and safely prepare proof of proper instrument decontamination. Confirmation of the process and the batch release step is done directly on the Careclave display without additional software. The protocols are saved on a USB memory stick or on the practice's network giving you the ability to fully verify and control the data at any time,
- Features and advantages of the MELAG Careclave 618 autoclave:
 - Lower service costs: Careclave is a much more affordable option than an autoclave and additional decontamination and care equipment. In addition, it eliminates costs associated with inspection and service of additional equipment,
 - 45% more space in the sterilization area: Save more space because there is no need for additional equipment,
 - Protection of instruments, reduction of servicing costs: Using the patented iDos dispensing process, each tool receives the right amount of oil. This ensures maximum reliability and tool life,
 - Revolutionary device concept: Not only does the Careclave enable the spectacular decontamination and care of up to 200 recessed dental instruments per day, but it also allows other instruments to be sterilized in the back of the chamber at the same time as the Carebox, in a single cycle. Indeed, Careclave meets the requirements of the standard for washer-disinfectors (PN EN ISO 15883) and autoclaves (PN EN 13060),
- Technical data:
 - Dimensions of the device (W x H x D): 48.0 x 50.0 x 65.3 cm height including the display 56.2 cm Careclave fits on a work surface with a depth of 60 cm,
 - Weight: 69 kg,
 - Chamber size / chamber volume: Ø 25 cm x 35 cm long / 18 liters,
 - Chamber loading quantity: max. 11 kg of instruments / 2 kg of textiles,
 - Carebox loading quantity: max 8 recessed instruments,
 - Batch sterilization: A and B,
 - Power supply: 220-240 V, 50 Hz,
 - Power consumption: 3 kW,
 - Compressed air connection: 4 - 8 bar (compressed air in accordance with ISO 7494-2),
 - Treated water supply: max. 35 °C, 1 - 10 bar, additional water treatment required,
 - Documentation interfaces: USB, network,
- Has a current Technical Passport issued,
- Warranty 12 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

