



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MELAG PREMIUM - PLUS Evolution Vacuklav 40 B+ Autoclave (New)**

- New,
- Made in Germany,
- The MELAG PREMIUM - PLUS Evolution Vacuklav 40 B+ autoclave is a model with a water supply connection,
- In addition to typical MELAG features such as time recording, intuitive operation, integrated documentation and load release, the steam sterilizers of the new Evolution series offer a wide range of unique and innovative solutions that could only have been refined to such perfection by a high-end medical hygiene and instrument cleaning specialist. The new Evolution series combines even shorter working times, wide instrument availability and energy savings
- Features and advantages of the MELAG PREMIUM - PLUS Evolution Vacuklav 40 B+ autoclave:
  - Time and energy savings: The patented DRYtelligence® technology is a unique and innovative solution from MELAG. This intelligent drying technology is based on a sensor-controlled drying function. A specially developed algorithm makes it possible to automatically adjust the drying according to the load in the sterilizer. DRYtelligence® technology makes it possible to automatically dry packaged and unpackaged instruments at full or incomplete charge, in a time period that depends on the fill level of the unit. DRYtelligence® not only guarantees optimal cleaning results, but also reduces cleaning time by up to 80%,
  - Work capacity: The batch volume of the Evolution series is an impressive 9 kg. This makes Evolution steam sterilizers among the most powerful in their category. In order to sterilize such a large load, we have equipped our sterilizers with new features to make the most of the chamber's capacity. Additional tray levels (up to eight trays) increase capacity, and also make it possible to sterilize individually packed handpieces and contra-angle handpieces, for example, in combination with MELAstore containers. In practice, this means fewer sterilization cycles, lower energy consumption and lower costs,
  - Secure documentation: the XXL color touchscreen display with a modern "flat design" enables intuitive operation at the highest level. The new Evolution series also features integrated software for documentation, release and tracking. Connecting the sterilizers to the MELAprint 60 label printer directly or via a computer enables the generation of customized documentation,
- Technical specifications:
  - Class B autoclave with fixed water supply connection,
  - Chamber size: D 35 cm x Ø 25 cm,
  - Chamber volume: 18 l,
  - Device dimensions (W x H x D): 46 x 50 x 55 cm,
  - Height with display: 56 cm,
  - Loading size: max. 9 kg of tools / 2 kg of textiles,
  - Requires a countertop with a minimum depth of 50 cm,
  - Weight: 55 kg,
  - Power supply: 220-230 V/50-60 Hz,
  - Power consumption: 3400 W,
  - Has a current Technical Passport issued,
- 12-month warranty,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

