



MELAG PREMIUM PLUS Evolution Vacuklav 41 B+ Autoclave (New)

- New,
- Made in Germany,
- Vacuklav 41B+ Evolution is a freestanding autoclave with an internal water tank,
- MELAG Vacuklav 41B+ Evolution is a powerful autoclave that has the same record-breaking fast sterilization cycles as the almost identical Vacuklav 40 B+ Evolution. Thanks to its innovative dual-chamber technology, it is capable of sterilizing packaged instruments in 20 minutes in a Class B program and unpackaged instruments in 10 minutes - including fractionated pre-vacuum and drying. DRYtelligence's patented intelligent drying shortens the cycle and saves resources. In addition, specially developed algorithms automatically adjust the drying time for each load. Thanks to its air-cooling capability, the Vacuklav 41 B+ Evolution is a true "free-standing" autoclave - only an electrical outlet is needed for installation. An important advantage of the MELAG Vacuklav 41B+ Evolution is its integrated software, which allows for release, labeling and tracking,
- Features and advantages of the MELAG PREMIUM - PLUS Evolution Vacuklav 41 B+ autoclave:
 - Fastest cycles with patented dual chamber technology,
 - Patented DRYtelligence system - intelligent drying for each load,
 - Documentation, release and tracking, large color touch screen with intuitive user interface,
 - Long sterilization chamber allowing batches with a maximum weight of 9 kg,
 - Water supply connections for cooling the powerful and reliable vacuum pump,
 - MELAconnect - device status and program progress via smartphone app,
- Technical specifications:
 - Class B freestanding autoclave with internal water tank,
 - Chamber size: G 35 cm x Ø 25 cm,
 - Chamber volume: 18 l,
 - Device dimensions (W x H x D): 46 x 50 x 55 cm,
 - Height with display: 56 cm,
 - Loading size: max. 9 kg of tools / 2 kg of textiles,
 - Requires a countertop with a minimum depth of 50 cm,
 - Weight: 60 kg,
 - Power supply: 220-230 V/50-60 Hz,
 - Power consumption: 3400 W,
- Has a current Technical Passport issued,
- Warranty 12 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

