



## MELAG PREMIUM PLUS Evolution Vacuklav 43 B+ Autoclave (New)

- New,
- Made in Germany,
- Autoclave MELAG PREMIUM - PLUS Evolution Vacuklav 43 B+ is a free-standing model with an internal water tank,
- MELAG Vacuklav 43B+ Evolution is an air-cooled autoclave is the ideal device for all clinics that need a fast, large and reliable autoclave without hydraulic connections. The deeper chamber (45 cm) allows sterilization of larger loads than in the similar Vacuklav 41 B+ Evolution. Very short sterilization times allow this autoclave to be used in clinics that want to sterilize numerous instruments in a short period of time. DRYtelligence's patented intelligent drying enables shorter working times and saves resources thanks to specially developed algorithms that adjust the drying time according to the load.
- Features and benefits of the MELAG PREMIUM - PLUS Evolution Vacuklav 43 B+ autoclave:
  - Fastest cycles with patented dual chamber technology,
  - Patented DRYtelligence system - intelligent drying for each load,
  - Documentation, release and tracking, large color touch screen with intuitive user interface,
  - Long sterilization chamber allowing batches with a maximum weight of 9 kg,
  - Free-standing autoclave with air cooling system,
  - MELAconnect - device status and program progress via smartphone app,
- Technical specifications:
  - Class B freestanding autoclave with internal water tank,
  - Chamber size: D 45 cm x Ø 25 cm,
  - Chamber volume: 23 l,
  - Dimensions of the device (W x H x D): 46 x 50 x 68 cm,
  - Height with display: 56 cm,
  - Loading size: max. 9 kg of tools / 2.5 kg of textiles,
  - Weight: 69 kg,
  - Power supply: 220-230 V/50-60 Hz,
  - Power consumption: 3400 W,
- Has a current Technical Passport issued,
- Warranty 12 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

