



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medical couch / diagnostic and treatment table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Power supply: 230 V,
- Upholstery with medical certificate (possibility to replace the upholstery with a new one from our color chart for an additional fee),
- The medical couch can also serve as a diagnostic and treatment table with the possibility of gynecological procedures,
- The large number of functions that the table provides allows for extended range of activities,
- Technical specifications:
  - Electric adjustment by foot control:
    - Height: 67 - 96 cm,
  - Manual adjustment:
    - Headrest backrest: -90° / 0° / 32°,
    - Back rest: 0° / 67°,
    - Seat: 0° / 20°,
  - Built-in side rails on both sides for attaching additional accessories,
  - Built-in paper rack,
  - Built-in 2 outlets designed to power additional accessories,
- Includes:
  - Foot controller,
  - Power cord,
  - Paper hanger,
- Dimensions:
  - Overall: 72 x 206 x 67/96 cm,
  - Headrest: 67 x 49 cm,
  - Shoulder rest: 67 x 35 cm,
  - Back rest: 67 x 82 cm,
  - Seat: 67 x 35 cm,
- Weight: 70 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



