



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LOVIBOND TC135S Thermostatic cabinet (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 270 VA,
- Technical specifications:
  - Temperature control in the range: from 2°C to 40°C, in steps of 0.1°C,
  - Capacity: 135 l,
  - Convenient internal 230V power outlets allow the operator to use mixing units or other electrical appliances without opening the door and lowering the temperature,
  - Built-in axial fans: 320 m3/h capacity,
  - Door with built-in lock,
- Included:
  - Physical instruction manual in GER, ENG, FR, ESP, PT, TR,
  - Power cord,
  - 3 pcs. - screen shelves,
  - Key,
- Dimensions:
  - Chamber: 51 x 44 x 70 cm,
  - Overall: 60 x 60 x 85 cm,
- Weight: 39 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

