



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LEISEGANG 2H Colposcope (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 40 VA,
- Bulb: Halogen 12 V / 20 W,
- Technical specifications:
  - Diopter adjustment: From - 6D to + 6D,
  - Free eyepiece spacing adjustment,
  - Infinitely adjustable light intensity,
  - Free side-to-side adjustment of the 2-segment arm and colposcope vertical tilt,
  - Wall mountable,
- Included:
  - Power cord,
- Dimensions: 134 x 71 cm,
- Weight: 11.5 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
  - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

