



## MEMMERT BE 700 Laboratory incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1800 W,
- Memmert BE700 incubator is electronically controlled and has a digital display with temperature control,
- Features:
  - PID microprocessor controller with selective digital setpoint specification for temperature and program times (on/off delay time),
  - Digital display of current temperature and various editing functions,
- Technical specifications:
  - Temperature control: from 30°C to 70°C,
  - Temperature fluctuation:  $\pm 0.5$  °C,
  - Integrated timer,
  - Capacity: 416 l,
  - Construction built of stainless steel, acid-resistant,
- Includes:
  - User manual in language (PDF),
  - Power cord,
- Dimensions:
  - Internal: 104 x 50 x 80 cm,
  - External: 119 x 65 x 106 cm,
- Weight: 125 KG,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

