



## **ALTHEA Nasopharyngoscope (New)**

- New,
- Realistic photos of the product,
- Made in Germany,
- Nasopharyngoscope ALTHEA is characterized by excellent, clear images and excellent ergonomics,
- The tip of the endoscope is tiltable and moves in two planes (within 180°), making the endoscope very handy during the procedure,
- Thanks to the use of modern technology, an unprecedented resolution of 18,000 pixels has been achieved,
- Technical specifications:
  - Working length: 300 mm,
  - Diameter: 2.8 mm,
  - Field of view: 75°,
  - Depth of field: 1 - 50 mm,
  - Tip bend: up: 150° / down: 150°,
- Included:
  - User manual,
  - Pump for leak test,
  - Original carrying case for transporting the instrumentation,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty: 24 months,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

