



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT BE 600 Laboratory incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1600 W,
- Technical specifications:
 - Infinitely adjustable temperature: from 20° to 70°C,
 - Infinitely adjustable time: from 1 minute to 999 hours or continuous operation mode,
 - Chamber interior made of stainless steel,
 - Chamber capacity: 256 l,
 - Chamber equipped with two sieve trays with possible arrangement,
 - Over-temperature alarm,
 - Built-in digital temperature display,
 - Adjustable fresh air intake,
 - Additional internal glass door,
 - Built-in RS 232 serial interface,
- Included:
 - Power cord,
- Dimensions:
 - Overall: 95 x 68.5 x 91.5 cm,
 - Internal: 80 x 50 x 64 cm,
- Weight: 87 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

