



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



**SAMSUNG UGEO G60 Ultrasound scanner (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Korea,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 620 VA,
- UGEO G60 Samsung Medison ultrasound scanner provides high working comfort and confident diagnosis with excellent quality diagnostic images,
- Technical specifications:
  - Types of examinations:
    - Gynecology: 3D / 4D,
    - Obstetrics: 3D / 4D,
  - Imaging modes:
    - 2D mode,
    - M mode,
    - Color DOPPLER,
    - Pulsed Wave (PW) Spectral Doppler,
    - Power Doppler (PD),
    - 3D/4D imaging modes,
    - Dual modes,
    - Quad modes,
    - Combined modes,
    - Simultaneous modes,
    - Zoom Mode,
  - Focusing:
    - Transmit focus, up to 8 points (4 points can be selected simultaneously),
    - Digital dynamic receive focus (continuous),
  - Display Technology: LCD,
  - Screen:
    - 18.5" in,
    - Resolution: 1366 x 768px,
  - Touchscreen:
    - 10.1" in,
    - Resolution: 1280 x 800px,
  - Probe connectors:
    - 4 active probe ports,
  - EVN4-9 (HS40) probe:
    - Frequency range: 4.0 - 9.0 MHz,
    - Application: obstetrics / gynecology / urology,
  - VE4-8 probe:
    - Frequency range: 4.0 - 8.0 MHz,
    - Application: abdomen / obstetrics / gynecology,
  - Inputs:
    - DVD / Ethernet / HDMI / USB 2.0 / VGA,
- Selectable language menu: EN, DE, FR, ES, IT, RU, CN, PT, BR,
- Included:
  - User manual in Polish,
  - Power cable,
  - BNC cable,
  - S-VIDEO cable,
  - Grounding cable,
  - EVN4-9 endovaginal probe,
  - Convex probe VE4-8,
- Overall dimensions: 55 x 98 x 160 cm,
- Weight: 120 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!  
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

