



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DATEX OHMEDA ADU S5 Anesthesia machine (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Sweden,
- Power supply: 230 V,
- Emergency power supply: through 4 new built-in 12 V / 2.3 Ah batteries,
- Datex Ohmeda ADU S5 uses electronic control of both gas feeding and vaporization, allowing full use of automatic data capture in managing care information and minimizing manual documentation,
- It has one port for Aladin cassettes along with one compatible vaporizer included,
- Specifications:
 - Ventilator measurement functions:
 - VT: 20 - 1400 mL,
 - f: 2 - 60 bpm,
 - Ti/Te: 1:4,5, 1:3, 1:2,5, 1:2, 1:1,5, 1:1, 2:1,
 - PEEP: 2 - 20 mbar,
 - The device has:
 - Compact patient breathing system with oxygen sensor, located on a movable arm,
 - Movable arm with attachment for patient tubing,
 - Removable DATEX OHMEDA Suction,
 - Rails for attaching additional accessories,
 - Tabletop with a working surface size of 44 x 30 cm,
 - Two gas cylinder holders,
 - Evac connector,
- Supplied with gases from the central installation: 2 x O₂, N₂O, 2 x AIR,
- Infinitely variable manual regulation of gas flow: O₂, AIR, N₂O (pressure gauges),
- Equipped with:
 - DATEX Cardiometer - OHMEDA Model S/5:
 - ComWheel navigation knob operation makes it easy to enter information and navigate the screen,
 - Custom display options allow the operator to view and analyze data in graphical and numerical formats, summarize measurements and trends in real time, and freely configure each element of the screen - from curves to parameter names,
 - Audible and visual alarms,
 - Screen diagonal: 10,5",
 - Overall dimensions: 33 x 26 x 16 cm,
 - DATEX OHMEDA S/5 ventilator monitor,
 - Inhalation agent vaporizer for DATEX OHMEDA ALADIN A-VSEVO sevoflurane:
 - Infusion type: Quik Fill,
 - Maximum capacity: 250 mL,
 - Modules:
 - M - ESTP: patient monitor module,
 - M - NIBP: pressure measurement module,
 - M - CAIOV: gas module,
- Mobile, on casters, including 2 with brake,
- Keypad on control panel in language: English,
- Included:
 - Gas supply hoses with tips of customer's choice (AGA or DIN),
 - User manual in English(PDF),
 - Power cord,
 - Patient breathing system,
 - Gas sampling line,
 - ECG 3-lead cable,
 - SPO₂ sensor,
- Dimensions: 100 x 75 x 141 cm,
- Weight: 163 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

