



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BTL 4825L Topline Laser and electrotherapy unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in U.K., (2012),
- Power supply: 230 V,
- Emergency power supply: through 2 built-in batteries 6 V / 1.2 Ah,
- Frequency: 50/60 Hz,
- Power: 72 W,
- BTL 4000 is used for electrotherapy treatments,
- Electrotherapy is another way of treating with electricity. It is a treatment during which various types of therapeutic currents are applied - direct or alternating, with different frequencies. Electrotherapy is a variation of physical therapy and is used in the symptomatic treatment of diseases of the musculoskeletal system, among others,
- Therapy parameters can be easily set using buttons on the touch screen and on the front panel of the device,
- Throughout the therapy, the device constantly displays on the screen information about the type of therapy used, accessories used, remaining therapy time and main parameters,
- The device has:
 - Ready-to-use therapy programs and user programs to help select the procedure,
 - A built-in encyclopedia with recommendations for the therapies used and many useful tips,
 - An extensive bank of diagnoses and programs,
 - Touch, easy-to-read 4.5" screen,
- Selectable menu in the following languages: EN, CS, DE, PL, ES, FR, IT, PT, RU, RO, TR, KR, NR, SE, HU, KZ, CN,
- Keypad on control panel in language: English,
- Included:
 - User manual in English (PDF),
 - Power cord and power supply,
 - IR infrared spot probe 400 mW / 830 nm,
 - Stylus,
 - 2 pcs. - BTL safety goggles,
 - 2 pcs. - patient cable,
 - 6 pcs. - BTL 70 x 50 mm electrode,
 - 2 pcs. - set of electrode cables,
 - 4 pcs. - electrode pouch,
 - 2 pcs. - viscose cover,
 - 2 pcs. - BTL electrodes,
 - Set of 4 pairs of attachment straps for BTL applicators,
- Dimensions: 34.5 x 16 x 14.5 cm,
- Weight: 4.5 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

